INT "G. PASCALE" NAPOLI

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CLINICAL CASE REVIEW 9

Bonopty®


Access for RF ablation of left scapula plasmacytoma

Case history

56 year old male with a history of multiple myeloma and intractable left shoulder pain, referred for RF ablation of known lytic lesions in the left scapula.

Procedure details

The 6.5 cm Bonopty® "short" Penetration Set (art no 10-1062) was used to access a superficial lytic lesion in the left scapula. The procedure was performed from a posterior approach with the patient in the prone position.

CT image showing multiple myeloma in left scapula.

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
Analysis of the procedure

Biopsy sample was not obtained at the time of the procedure. The lesion had been previously sampled. Once access was gained with the 6.5 cm Bonopty® "short" Penetration Set, the therapeutic RF ablation procedure was performed without difficulty.

Comments

The shorter Bonopty® 6.5 cm Penetration Set (art no 10-1062) provided proper access to the bone and the stability needed for the ablation procedure of the superficial lytic lesion in the left scapula.

Bonopty® "short" Penetration Set (art no 10-1062) accessing the lesion and creating a path for ablation probe.

2007 Rev. 01/2011 it

Head Office:
AprioMed AB
Viridings Alle 28
SE-751 50 Uppsala, Sweden
Phone: +46 18 430 14 40
Free Fax (Int'l): +46 770 339972

E-mail: info@apriomed.com
www.apriomed.com

US Sales Office:
AprioMed, Inc.
2 Palmer Dr, Suite 1a
Londonderry, NH 03053, USA
Phone: +1 603 421 0875
Fax (US toll free): 866 239 1797

CLINICAL CASE REVIEW 1

Bonopty®

Penetration of 8 mm cortical bone

Case history

Seventy-five year old man with a suspected metastasis in the femoral shaft.

Biopsy details


Penetration of the 8mm thick, cortical bone was achieved by manually drilling with the Bonopty® Drill. The drill was then exchanged for the Bonopty® Biopsy Cannula which was used to obtain two samples in the bone marrow, through the anchored Bonopty® Penetration Cannula.

Analysis of the samples

Histologic examination of these samples revealed metastatic bone marrow involvement.

Comments

The intact, thick cortical bone was easily penetrated with the eccentric drill using only finger-tip force.

CT scan of the femur

The anchored Bonopty® Penetration Cannula

Samples obtained with the Bonopty® Biopsy Cannula

8 mm thick cortical bone

Approved by: Dr. K. Gustav, U. Åstrand, Dr. S. Hult, M. Åberg, Department of Orthopedic Hospital, Uppsala University Hospital, Uppsala, Sweden.

Copyright © 2000

Head Office:
AprioMed AB
Vindings Allé 28
SE-754 56 Uppsala, Sweden

Phone: +46 18 430 14 10
Free Fax (Int'l): +46 770 339972
E-mail: info@apriomed.com
www.apriomed.com

AprioMed

Gangi, A.; Guth, S.; Diemann, J.; Roy, C. (2001). Interventional musculoskeletal procedures. *RadioGraphics*, 21(2), e1.

Kjar, R. A.; Powell, G. J.; Schlicht, S. M.; Smith, P. J.; Slavin, J.; Choong, P. F. (2006). Percutaneous radiofrequency ablation for osteoid osteoma: experience with a new treatment. *Medical Journal of Australia*, 184(11), 563-565.

OBJECTIVE: To review our experience with percutaneous radiofrequency ablation (RFA) for osteoid osteoma.

PARTICIPANTS, DESIGN AND SETTING: Retrospective review of 24 patients with osteoid osteoma treated with percutaneous RFA at St Vincent's Hospital, Melbourne, from August 2000 to February 2005.

MAIN OUTCOME MEASURES: Initial response to treatment, return of symptoms, time to recurrence, complications and histopathological correlation. **RESULTS:** In 23 of 24 patients, there was immediate relief of symptoms. One-third of patients had a return of symptoms. The 24 patients underwent a total of 32 RFA procedures. Of the eight patients who had a recurrence, five had an initial lesion ≥ 10 mm in maximum diameter. Twenty-three of the 24 patients were pain-free at the end of the study period. Patients were followed up for a median of 26 months. There were no long-term complications.

CONCLUSIONS: Percutaneous RFA is a safe and efficacious treatment for osteoid osteoma with a low morbidity rate. Despite recurrence after treatment, re-treatment is simple and effective.

Morrison, W.; Sanders, T. (2008). *Problem solving in musculoskeletal imaging*. Elsevier.

Peterson, J.; Fenton, D.; Czervionke, L. (2007). *Image-guided musculoskeletal intervention*. Elsevier.

Pinto, C. H.; Taminiau, A. H.; Vanderschueren, G. M.; Hogendoorn, P. C.; Obermann, W. R. (2002). Technical considerations in CT-guided radiofrequency thermal ablation of osteoid osteoma: tricks of the trade. *American Journal of Roentgenology*, 179(6), 1633-1642.

Roberts, C. C.; Morrison, W. B.; Leslie, K. O.; Carrino, J. A.; Lozevski, J. L.; Liu, P. T. (2005). Assessment of bone biopsy needles for sample size, specimen quality and ease of use. *Skeletal Radiology*, 34(6), 329-335.


Rosenthal, D.I.; Hornicek, F. J.; Torriani, M.; Gebhardt, M. C.; Mankin, H. J. (2003). Osteoid osteoma: percutaneous treatment with radiofrequency energy. *Radiology*, 229(1), 171-175.

Rosenthal, D.I.; Marota, J. J.; Hornicek, F. J. (2003). Osteoid osteoma: elevation of cardiac and respiratory rates at biopsy needle entry into tumor in 10 patients. *Radiology*, 226(1), 125-128.

Tan, K.; Yoong, P.; Marshall, T. J.; Martin, C. (2012). Percutaneous drainage as a novel approach for the treatment of Brodie's abscess. *Clinical Radiology*, 67(10), 1030-1033.

Vanderschueren, G. M.; Taminiau, A. H.; Obermann, W. R.; van den Berg-Huysmans, A. A.; Bloem, J. L. (2004). Osteoid osteoma: factors for increased risk of unsuccessful thermal coagulation. *Radiology*, 233(3), 757-762.

Yip, P. S.; Lam, Y. L.; Chan, M. K.; Shu, J. S.; Lai, K. C.; So, Y. C. (2006). Computed tomography-guided percutaneous radiofrequency ablation of osteoid osteoma: local experience. *Hong Kong Medical Journal*, 12(4), 305-309.

Bonoptoy® Adapter 1:3 Revolutions

Adapter 1:3 Revolutions

Universal connector, fits both 12G and 14G drills

Sterile, single use

Ordering information:

250 Bonoptoy® Adapter 1:3 Revolutions

MOQ 5 units

DOCUMENTS

Product Sheets (PDF) ▾

Presentations (PDF) ▾

Movies (YouTube) ▾

Case studies (PDF) ▾

Articles and references ▾

Click reference to view abstract

Ahlstrom, K. H.; Astrom, K. G. (1993). CT-guided bone biopsy performed by means of a coaxial biopsy system with an eccentric drill. *Radiology*, 188(2), 549-552.

Astrom, K. G.; Ahlstrom, K. H.; Hagberg, H. E. (1997). CT-guided transvertebral core biopsy of a retrocrural mass. *American Journal of Roentgenology*, 169(4), 991-993

Astrom, K. G.; Sundstrom, J. C.; Lindgren, P. G.; Ahlstrom, K. H. (1995). Automatic biopsy instruments used through a coaxial bone biopsy system with an eccentric drill tip. *Acta Radiologica*, 36(3), 237-242.

Ciray, I.; Astrom, G.; Sundstrom, C.; Hagberg, H.; Ahlstrom, H. (1997). Assessment of suspected bone metastases. CT with and without clinical information compared to CT-guided bone biopsy. *Acta Radiologica*, 38(5), 890-895.

BONOPTY® 14G

Coaxial bone biopsy system

The Bonopty® 14G system for bone biopsies has been a success since it first reached the market. Today, Bonopty® 14G exists in three lengths, all to perfectly match your various needs.

- Gain access into the bone, even through thick cortical bone
- Coaxial system facilitates multiple sampling, or treatment of lesions
- Achieve excellent core samples with few crush artifacts
- Successfully sample somewhat sclerotic lesions
- A variety of lengths enhance perfect match for every procedure

PRODUCT ASSORTMENT

Bonopty® Penetration Set, 14G

Bonopty® Penetration Set, 14G

Penetration Cannula with Stylet (ID 1.8 mm, OD 2.1 mm, Length 6.5/9.5/13 cm)

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ADVANCED MEDICAL SUPPLIES S.p.A.